



SPA463Hu02 100µg
Synthetic Thyroid Stimulating Hormone (TSH)
Organism Species: *Homo sapiens* (Human)
Instruction manual

FOR RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

[**PROPERTIES**]

Source: Synthetic Peptide

Residues: Glu58~Tyr74

Purity: >90%

Traits: Freeze-dried powder

Buffer Formulation: 10mM PBS, pH7.4.

Predicted isoelectric point: 6.2

Predicted Molecular Mass: 2180Da

Applications: Immunogen; Blocking Peptide.

[**USAGE**]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[**STORAGE AND STABILITY**]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[**SEQUENCE**]

ERR ECAYCLTINT TICAGY

[IDENTIFICATION]




Figure 1. HPLC




Figure 2. MASS Spectrometry

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.