

Featured Products From Cloud-Clone Corp. — Immunity to infection

Name	Protein	Antibody	Assay kit
High Mobility Group Protein 1 (HMG1)	RPA399Ca01, RPA399Hu01, RPA399Mu01, RPA399Po01, RPA399Ra01, RPA399Mu02, RPA399Hu02, RPA399Hu61 ...	PAA399Ca01, PAA399Hu01, MAA399Hu22, PAA399Mu01, PAA399Po01, PAA399Ra01, PAA399Hu71, PAA399Mu71, PAA399Ra71, PAA399Hu81, PAA399Mu81, PAA399Ra81, PAA399Mu02, MAA399Ca21, MAA399Po21, MAA399Ra21, PAA399Hu02, PAA399Hu06. ...	SEA399Hu, SEA399Ra, SEA399Mu, SEA399Ca, SEA399Po, SCA399Hu ...
C Reactive Protein (CRP)	RPA821Ra01, RPA821Rb01, RPA821Po01, RPA821Gu01, RPA821Mu01, RPA821Bo01, RPA821Ga01, RPA821Mu02, RPA821Rb02, RPA821Hu01, NPA821Hu02. ...	PAA821Mu01, PAA821Ra01, MAA821Hu22, MAA821Mu22, MAA821Ra22, PAA821Gu01, PAA821Rb52, PAA821Po01, PAA821Bo01, PAA821Ga01, PAA821Fi01, PAA821Rb51, MAA821Po22, MAA821Ra82, MAA821Po82, MAA821Ra72, MAA821Po72, PAA821Gu71, PAA821Hu71, PAA821Mu71, PAA821Rb71, PAA821Ra71, PAA821Po71, PAA821Gu81, PAA821Mu81, PAA821Rb81, PAA821Ra81, PAA821Po81, PAA821Mu02, MAA821Bo21, MAA821Ga21, MAA821Gu21, PAA821Hu01, PAA821Hu02, MAA821Hu72. ...	SEA821Hu, SEA821Bo, SEA821Ra, SEA821Mu, SEA821Rb, SEA821Ga, SEA821Po, SEA821Gu, SCA821Hu, SCA821Ra, AEA821Hu, HEA821Hu, HEA821Ra. ...
Complement Component 3a (C3a)	RPA387Mu01, RPA387Ra01, RPA387Hu01, RPA387Po01 ...	PAA387Mu01, PAA387Ra01, PAA387Hu01, PAA387Po01. ...	SEA387Hu, SEA387Ra, SEA387Mu, SEA387Ca, SCA387Hu, SEA387Po, SEA387Rb. ...
Complement Component 5a (C5a)	RPA388Mu01, RPA388Ra01, RPA388Ra02, RPA388Hu01, RPA388Hu61. ...	PAA388Mu01, PAA388Ra01, PAA388Ra02, PAA388Ra81, PAA388Ra71, PAA388Mu71, PAA388Hu01, PAA388Hu06. ...	SEA388Hu, SEA388Ra, SEA388Mu, SEA388Po, SEA388Rb, SEA388Ca, SCA388Mu, SCA388Ra. ...
Tryptase (TPS)	RPB070Ca01, RPB070Hu01, RPB070Mu01, RPB070Po01, RPB070Ra01, RPB070Ra02, RPB070Eq01, RPB070Mu02. ...	PAB070Ca01, PAB070Hu01, MAB070Hu22, PAB070Mu01, PAB070Po01, PAB070Ra01, PAB070Ra02, PAB070Ca71, PAB070Hu71, PAB070Mu71, PAB070Ra71, PAB070Eq01, PAB070Mu02, MAB070Ca21, MAB070Eq21, MAB070Po21. ...	SEB070Hu, SEB070Po, SEB070Ra, SEB070Ca, SEB070Ra, SEB070Mu, SEB070Eq, SCB070Hu. ...
Defensin Alpha 1, Neutrophil (DEFA1)	RPB705Hu01, RPB705Mu01, RPB705Mu02. ...	PAB705Hu01, MAB705Hu22, PAB705Mu01, PAB705Mu02. ...	SEB705Hu, SEB705Mu ...
Procalcitonin (PCT)	RPA689Mu01, RPA689Hu01, RPA689Ga01, RPA689Ra01 ...	PAA689Hu01, PAA689Mu01, PAA689Ra01, MAA689Hu21, MAA689Ra21, PAA689Hu71, PAA689Ga01. ...	SEA689Hu, SCA689Hu, SEA689Ra, SEA689Mu ...
Lipopolysaccharide Binding Protein (LBP)	RPB406Bo01, RPB406Hu01, RPB406Mu01, RPB406Ra01, RPB406Mu02, RPB406Po01, RPB406Rb01. ...	PAB406Bo01, PAB406Hu01, MAB406Hu22, PAB406Mu01, PAB406Ra01, PAB406Hu71, PAB406Mu02, PAB406Po01, MAB406Bo21, MAB406Po21, MAB406Ra21, PAB406Rb51, PAB406Bo71. ...	SEB406Hu, SEB406Ra, SEB406Mu, SEB406Bo, SEB406Po. ...
Toll Like Receptor 2 (TLR2)	RPA663Hu01, RPA663Mu01, RPA663Ra01, RPA663Hu02, RPA663Mu02. ...	PAA663Hu01, PAA663Mu01, PAA663Ra01, PAA663Hu02, PAA663Hu71, PAA663Hu81, MAA663Hu21, PAA663Mu02, MAA663Ra21, MAA663Hu22. ...	SEA663Ra, SEA663Hu, SEA663Mu. ...
Alpha-1-Antitrypsin (a1AT)	RPB697Hu01, RPB697Mu01, RPB697Ra01, RPB697Mu02 ...	PAB697Hu01, PAB697Mu01, PAB697Ra01, PAB697Rb51, PAB697Po01, MAB697Hu22, PAB697Hu71, PAB697Mu71, PAB697Ra71, PAB697Mu81, PAB697Mu02, MAB697Ra21. ...	SEB697Hu, SEB697Mu, SEB697Ra, SCB697Hu ...
Transglutaminase 1, Keratinocyte (TGM1)	RPB773Hu01, RPB773Rb01, RPB773Ra01, RPB773Ca01, RPB773Mu01, RPB773Mu02. ...	PAB773Hu01, MAB773Hu22, PAB773Ra01, PAB773Ca01, PAB773Mu01, PAB773Ca81, PAB773Hu81, PAB773Mu81, PAB773Ca71, PAB773Hu71, PAB773Mu71, PAB773Rb51, PAB773Mu02. ...	SEB773Hu, SEB773Mu ...
S100 Calcium Binding Protein A8 (S100A8)	RPB792Bo01, RPB792Hu01, RPB792Mu01, RPB792Ra01, RPB792Mu02. ...	PAB792Bo01, PAB792Hu01, PAB792Mu01, PAB792Ra01, PAB792Bo71, PAB792Hu71, PAB792Mu71, PAB792Ra71, PAB792Bo81, PAB792Hu81, PAB792Mu81, PAB792Ra81, MAB792Hu21, PAB792Mu02, MAB792Bo21, MAB792Ra21. ...	SEB792Hu, SEB792Mu, SEB792Bo, SEB792Ra ...
S100 Calcium Binding Protein A9 (S100A9)	RPB793Bo01, RPB793Hu01, RPB793Mu01, RPB793Ra01, RPB793Mu02. ...	PAB793Bo01, PAB793Hu01, MAB793Hu22, PAB793Mu01, PAB793Ra01, MAB793Hu82, MAB793Hu72, PAB793Hu71, PAB793Mu71, PAB793Ra71, PAB793Hu81, PAB793Mu81, PAB793Mu02, MAB793Bo21. ...	SEB793Hu, SEB793Mu, SEB793Bo. ...
Alpha-2-Heremans Schmid Glycoprotein (aHSG)	RPA178Bo01, RPA178Cp01, RPA178Gu01, RPA178Hu01, RPA178Mu01, RPA178Po01, RPA178Ra01, RPA178Rb01, RPA178Si01, RPA178Hu02, RPA178Mu02. ...	PAA178Hu01, PAA178Mu01, PAA178Ra01, PAA178Gu01, PAA178Rb51, PAA178Si01, PAA178Po01, PAA178Bo01, PAA178Cp01, MAA178Hu22, PAA178Hu71, PAA178Hu02, PAA178Mu02, MAA178Bo21, MAA178Po21, MAA178Ra21. ...	SEA178Hu, SEA178Mu, SEA178Po, SEA178Ra, SEA178Bo. ...
Gelsolin (GS)	RPA372Hu01, RPA372Mu01, RPA372Ra01, RPA372Si01, RPA372Bo01, RPA372Mu02, RPA372Hu02, RPA372Hu03. ...	PAA372Hu01, PAA372Mu01, PAA372Ra01, PAA372Si01, MAA372Hu22, PAA372Hu71, PAA372Ra71, PAA372Ra81, PAA372Bo01, PAA372Mu02, PAA372Hu02, PAA372Hu03. ...	SEA372Hu, SEA372Mu, SCA372Hu. ...
Cathelicidin Antimicrobial Peptide (CAMP)	RPC419Hu01, RPC419Mu01, RPC419Ra01, RPC419Mu02, CPC419Hu21. ...	PAC419Hu01, MAC419Hu22, PAC419Hu71, PAC419Mu01, PAC419Ra01, PAC419Mu02, PAC419Hu02. ...	SEC419Mu, CEC419Hu ...

For more products, please check <http://www.cloud-clone.us/>.



Cloud-Clone Corp.



Name	Protein	Antibody	Assay kit
Serum Amyloid A (SAA)	RPA885Bo01,RPA885Hu01, RPA885Mu01,RPA885Rb01, RPA885Mu02,RPA885Rb02, CPA885Hu21...	PAA885Bo01,PAA885Hu01,MAA885Hu22,PAA885Mu01, PAA885Rb52,PAA885Rb51,PAA885Hu71,PAA885Mu71, PAA885Rb71,PAA885Hu81,PAA885Mu81,PAA885Rb81, PAA885Mu02,MAA885Bo21,PAA885Hu02...	SEA885Mu,SEA885Bo, SEA885Rb,SEA885Hu, SCA885Rb,CEA885Hu ...
Clara Cell Protein 16 (CC16)	RPA857Hu01,RPA857Mu01, RPA857Ra01,RPA857Mu02 ...	PAA857Hu01,MAA857Hu22,PAA857Mu01,PAA857Ra01, PAA857Hu71,PAA857Mu71,PAA857Ra71,PAA857Hu81, PAA857Mu81,PAA857Ra81,PAA857Mu02,MAA857Ra21...	SEA857Hu,SEA857Ra, SEA857Mu...
Aspartate Aminotransferase (AST)	RPB214Hu01,RPB214Mu01, RPB214Ra01,RPB214Bo01, RPB214Mu02,RPB214Hu51 ...	PAB214Hu01,MAB214Hu22,PAB214Mu01,PAB214Ra01, PAB214Hu71,PAB214Mu71,PAB214Ra71,PAB214Hu81, PAB214Mu81,PAB214Bo01,PAB214Mu02,MAB214Bo21, MAB214Ra21,PAB214Hu51,PAB214Hu07...	SEB214Ra,SEB214Mu, SEB214Hu,SEB214Bo ...
Alanine Aminotransferase (ALT)	RPA207Bo01,RPA207Gu01, RPA207Hu01,RPA207Mu01, RPA207Po01,RPA207Ra01, RPA207Rb01,RPA207Hu02, RPA207Hu03,RPA207Ra02, RPA207Mu02...	PAA207Hu01,PAA207Mu01,PAA207Ra01,PAA207Bo01, PAA207Gu01,MAA207Hu22,PAA207Po01,PAA207Rb51, PAA207Hu02,PAA207Po81,PAA207Po71,PAA207Hu03, PAA207Hu71,PAA207Ra71,PAA207Hu81,PAA207Ra81, PAA207Ra02,MAA207Ra21,PAA207Mu02,MAA207Bo21, MAA207Hu72...	SEA207Mu,SEA207Hu, SEA207Ra,SCA207Ra, SEA207Bo...
Elastase 2, Neutrophil (ELA2)	RPA181Ca01,RPA181Hu01, RPA181Mu01,RPA181Ra01, RPA181Mu02...	PAA181Ca01,PAA181Hu01,MAA181Hu22,PAA181Mu01, PAA181Ra01,PAA181Mu02,MAA181Ca21,MAA181Ra21 ...	SEA181Ra,SEA181Hu, SEA181Mu,SEA181Ca ...
Defensin Alpha 5, Paneth Cell Specific (DEFA5)	RPB912Hu01,RPB912Mu01, RPB912Mu02...	PAB912Hu01,MAB912Hu22,PAB912Hu71,PAB912Mu01, PAB912Mu02...	SEB912Mu,SEB912Hu ...
Macrophage Migration Inhibitory Factor (MIF)	RPA698Hu01,RPA698Mu01, RPA698Ra01,RPA698Mu02, RPA698Hu02...	PAA698Hu01,MAA698Hu22,PAA698Mu01,PAA698Ra01, PAA698Mu02,MAA698Ra21,PAA698Hu02...	SEA698Hu,SEA698Mu, SEA698Ra,SCA698Mu ...
Immunoglobulin G (IgG)	NPA544Hu01,NPA544Mu01, NPA544Bo01,NPA544Rb01, NPA544Ca01,NPA544Cp01, NPA544Eq01,NPA544Ga01, NPA544Gu01,NPA544Ov01, NPA544Po01,NPA544Ra01, RPA544Rb61...	PAA544Hu01,PAA544Mu01,PAA544Ra01,PAA544Bo01, PAA544Po01,PAA544Ga01,PAA544Cp01,PAA544Ca01, PAA544Gu01,PAA544Rb52,PAA544Fe01,PAA544Rb51, PAA544Si01,MAA544Bo22,PAA544Ra02,MAA544Ra22, MAA544Ov22,MAA544Eq22,MAA544Ov82,PAA544Ca91, MAA544Eq82,MAA544Ov72,MAA544Eq72,PAA544Bo71, PAA544Ca71,PAA544Ga71,PAA544Fe71,PAA544Hu71, PAA544Mu71,PAA544Rb71,PAA544Si71,PAA544Po71, PAA544Gu81,PAA544Ga81,PAA544Mu81,PAA544Ra81, MAA544Ra21,MAA544Hu22,MAA544Gu21, MAA544Bo21,MAA544Rb21...	CEA544Gu,CEA544Ga, CEA544Po,CEA544Rb, CEA544Hu,CEA544Mu, CEA544Ra,CEA544Ca, CEA544Bo,CEA544Co, CCA544Ca,CCA544Cp, CCA544Eq,CCA544Ga, CCA544Gu,CCA544Hu, CCA544Mu,CCA544Po, CCA544Ra,CCA544Rb, CEA544Cp,CEA544Eq, HEA544Bo,HEA544Hu, HEA544Mu,SEA544Cp...
Surfactant Associated Protein D (SPD)	RPB039Hu01,RPB039Mu01, RPB039Ra01,RPB039Mu02 ...	PAB039Hu01,MAB039Hu22,PAB039Mu01,PAB039Ra01, PAB039Mu02,PAB039Hu71,PAB039Mu71,PAB039Ra71, PAB039Hu81,PAB039Ra81,MAB039Ra21...	SEB039Hu,SEB039Ra, SEB039Mu...
Alpha-2-Macroglobulin (a2M)	RPB017Hu01,RPB017Mu01, RPB017Ra01,RPB017Ra02, RPB017Ra03,RPB017Mu02, RPB017Mu03,RPB017Po01, RPB017Rb01...	PAB017Hu01,PAB017Mu01,PAB017Ra01,PAB017Rb51, PAB017Po01,MAB017Hu22,PAB017An01,PAB017Ra02, PAB017Mu02,PAB017Mu03,PAB017Ra03,PAB017An71, PAB017Hu71,PAB017Ra71,MAB017Ra21...	SEB017Hu,CEB017Mu, SCB017Hu,HEB017Hu, SEB017Ra...
Secretory Leukocyte Peptidase Inhibitor (SLPI)	RPB312Hu01,RPB312Mu01, RPB312Mu02...	PAB312Hu01,MAB312Hu22,PAB312Mu01,PAB312Mu71, PAB312Mu81,PAB312Mu02...	SEB312Hu,SEB312Mu, SCB312Mu...
Regenerating Islet Derived Protein 3 Gamma (REG3g)	RPE676Hu01,RPE676Mu01, RPE676Ra01,RPE676Mu02 ...	PAE676Hu01,PAE676Mu01,PAE676Mu71,PAE676Mu81, PAE676Ra01,MAE676Hu21,PAE676Mu02...	SEE676Hu,SEE676Mu, ...
Leukocyte Elastase Inhibitor (LEI)	RPD387Mu01,RPD387Hu01, RPD387Mu02,RPD387Ra01...	PAD387Mu01,PAD387Hu01,PAD387Mu02,MAD387Hu21, PAD387Ra01...	SED387Mu,SED387Hu ...
Fibronectin (FN)	RPA037Hu01,RPA037Mu01, RPA037Ra01,RPA037Ga01, RPA037Hu02,RPA037Bo01, RPA037Mu02,RPA037Po01, NPA037Bo01,RPA037Hu03 ...	PAA037Hu01,PAA037Mu01,PAA037Ra01,PAA037Gu01, PAA037Ca01,PAA037Po01,PAA037Bo01,PAA037Cp01, PAA037Eq01,MAA037Hu22,PAA037Ga01,PAA037Hu02, PAA037Bo71,PAA037Ga71,PAA037Hu71,PAA037Mu71, PAA037Ra71,PAA037Ga81,PAA037Ra81,MAA037Bo21, MAA037Po21,MAA037Ra21,PAA037Mu02,PAA037Hu03, PAA037Hu04...	SEA037Hu,SEA037Mu, SEA037Ra,SEA037Bo, SCA037Hu,SEA037Bo, SEA037Po,HEA037Hu ...
Vanin 1 (VNN1)	RPC598Hu01,RPC598Ra01, RPC598Mu01,RPC598Hu02, RPC598Hu03,RPC598Bo01...	PAC598Hu01,MAC598Hu22,PAC598Ra01,PAC598Mu01, PAC598Hu02,PAC598Hu03,PAC598Hu71,PAC598Mu71, MAC598Ra21,PAC598Bo01...	SEC598Ra,SEC598Hu ...
Arginase (Arg)	RPB120Hu01,RPB120Mu01, RPB120Ra01,RPB120Mu02 ...	PAB120Hu01,PAB120Mu01,PAB120Ra01,MAB120Hu22, PAB120Rb51,PAB120Hu71,PAB120Mu71,PAB120Rb71, PAB120Ra71,PAB120Hu81,PAB120Mu81,PAB120Ra81, MAB120Hu21,PAB120Mu02,MAB120Ra21,MAB120Hu71 ...	SEB120Hu,SEB120Mu, SEB120Ra...

For more products, please check <http://www.cloud-clone.us/>.